medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennnes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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medi Trading (Shanghai) Co. Ltd.
Room 613, HnR Group Building
No.1178, Tian Yao Qiao Road
Xuhui District, Shanghai, 200232
People's Republic of China
T +86-21 50582318
info@medi.cn
www.medi.cn

Medi Turk Ortopedi Medikal İc ve Dis Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 16 /4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi Ukraine LLC
Evhen Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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circaid® juxtacures®

Intended purpose

This compression product serves to compress  the leg in patients with venous and lymphatic diseases.

Indications

- Chronic venous insufficiency
 - Varicose Veins
 - Lipodermatosclerosis
 - Ulcus cruris venosum
 - Mixed arterial-venous ulcer (Ulcus cruris mixtum)
- Deep vein thrombosis/Thrombosis Prevention
- Acute DVT
- Post Thrombotic Syndrome
- Post Sclerotherapy
- Dependent oedema

Contraindications

- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5; ankle arterial pressure < 60mmHg; toe pressure < 30 mmHg or TcPO2 < 20 mmHg on dorsum of foot). If non-elastic materials are used, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot
- Any circumstances under which an increase in venous or lymphatic reflux is undesirable
- Suspected or known untreated acute deep vein thrombosis

Risks / Side effects

- Slight or moderate peripheral arterial disease
- Impaired sensation – patient must be

- able to feel the degree of compression
- Sensitivity to one of the used materials (may cause skin irritation or itching)
- Paediatric and assisted use – patients must be able to sense and communicate applied pressure

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Note on wearing

When applied correctly, the straps overlap slightly so that there are no gaps. Only tighten the straps in the intended places. The hook-and-loop fastener must never come into contact with the skin. The compression product should feel firm, but comfortable. Unless your treating physician orders or recommends otherwise, we recommend wearing the legwear together with the respective footwear included in the set. If in doubt, please consult your treating physician.

Individualizing

- Separate the body and spine from one another, unroll the bottom band and locate the number markings along the bottom and top edge of the body piece.
- Locate the patient's ankle circumference measurement amongst the printed numbers along the bottom edge. Then, locate the patient's calf circumference measurement amongst the printed numbers along the top edge of the garment.
- Take the spine and place the bottom hook tab on the appropriate number (2a). Take the top hook tab of the spine and place it on the appropriate number (2b). Adjust the center spine tab to create a smooth spine curve (2c).

Note: *There is no one correct spine shape. The garment may not lay flat once the spine has been set (2d). This is to be expected when dealing with more extreme spine angles. If there is a large difference between the two circumferences, the spine may not reach both the bottom and top edges of the body piece. In this case,*

center the spine (2e). When dealing with more extreme spine angles, excess material sticking out beyond the top of the spine piece can be trimmed (2f).

- With the spine in place, fold over the smaller garment piece. Cut the excess material from the body piece leaving approximately 5 cm of overlap (3). Be sure not to cut any other part of the garment.

Note: *If desired, attach the separate hook tab stays along the spine for additional adhesion (4).*

Tip: *You may also trace the spine to easily reattach it in its original position in the unlikely event of the spine becoming detached.*

- Take the BPS-card and highlight the range that contains the patient's ankle circumference. If the patient's measurement falls between two ranges on the card, 30.5 cm for example, choose the smaller range.

Fitting instructions

- Slide the circaid undersleeve onto the lower leg and over the knee. The cuffed edge of the undersleeve should cover the heel. Ensure that there are no wrinkles in the undersleeve (5a).
- Align the garment just above the ankle bone and below the knee crease. The black side of the material should face the skin (5b).
Unroll the top two bands and loosely secure them to the outside of the garment to anchor it in place (5c). This enables you to easily raise or lower the garment for correct positioning.

Note: *The bottom most band can be identified by the circaid-logo on the beige side of the material.*

- Unroll the bottom two bands and secure the bottom most band to the outside of the legging to a firm and comfortable compression level (5d).
- While holding the second band, detach the next band (5e). Secure the second band and continue this process with the other bands, going up the leg (5f).
- The garment should lay flat and wrinkle-free against the leg. Adjust each band to a firm and comfortable compression level (5g). Ensure that the bands are placed in an alternating order.
- Once all bands are secured, inspect the garment for any gaps or creases. Adjust the bands as necessary.

Using the built-in-pressure system (BPS)

The built-in-pressure system (BPS) is im-

plemented on the circaid product with the BPS card and the BPS lines and is used to set, check and adjust the defined compression range of the circaid compression product. You can use the BPS card to set the prescribed compression range.

- Put on the circaid compression product so that it fits tightly, but comfortably (see fitting instructions).
- Pick up the BPS card.
- Locate your ankle circumference on the card and mark the corresponding edge for your reference. In addition, identify the degree of compression prescribed for you on the colour coded card.
- Start with the lowest strap on the compression product. Place the black triangle on the BPS card against one of the BPS lines; the second BPS line is positioned on the prescribed degree of compression.
- If the BPS line is below (B) or above (C) the desired compression area, loosen and fasten the strap again until reaching the desired compression (A).
- Repeat placing the BPS card and setting the compression range for all straps from distal to proximal. (5h)
Always make sure that all straps are set to the same, prescribed level of compression to ensure graduated compression.

In order to maintain the prescribed compressive pressure throughout the day, check the distances between the BPS lines with the BPS card and adjust the straps at intervals throughout the day, if necessary. Always make sure that the correct, prescribed pressure is applied for compression.

Note: The BPS lines for individual straps are generally located at different points; it is not necessary to align the BPS lines in one row lengthwise or in any other specific position.

Note: Once the correct pressure is set (5h), you may fold the top and bottom part of the circaid undersleeve over the top and bottom edge of the circaid juxtacures garment (5i).

Reduction of the leg circumference

Remove the patient's garment and remeasure the ankle and calf circumference. Adjust the position of the small part of the garment accordingly (see Individualizing).

Once the device has been properly adjusted, re-apply the legging to the patient's limb. Be sure to also check whether or not the patient should be using a different edge of the BPS-card as the ap-

propriate range on the card may have changed. Make the appropriate adjustments, if necessary.

Tip Once the patient's ulcer has healed, circaid juxtalite and mediven round-knitted stockings are available for order to prevent ulcer recurrence.

Care instructions

This compressive product can be washed on a delicate cycle and dried in a tumble dryer set to a low heat. Please roll up the straps for laundry, the hook-and-loop-fasteners must not be open. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

-
-  Machine wash warm
 -  Do not bleach
 -  Tumble dry, low heat
 -  Do not iron
 -  Do not dry clean
 -  Do not wring
-

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Lifespan

6 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product). Unless otherwise prescribed by your doctor, wear the appliance daily from morning to night.

Material composition

Nylon, elastane, PU

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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