medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
www.medi.de

medi Australia Pty Ltd
Unit 4/13-21 Hallmark Street
Pendle Hill NSW 2145
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@medi.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
104-1375 Lionel-Boulet,
Varennnes, Québec,
QC Canada J3X 1P7
T +1 450-583-3317 / +1 800 361 3153
F +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

UK Responsible Person
medi UK Ltd.
Plough Lane
Hereford HR4 0EL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi

circaid® juxtafit® essentials / premium lower leg

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Brugsvejledning. Bruksanvisning. Návod k použití. Инструкция по использованию. Instrukcja zakładania. Használati útmutató. Інструкції з використання. دليل الاستخدام. Brugsanvisning. Kasutusjuhend. הוראות שימוש. Käyttöohje. 使用説明書.

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


Important notes

This medical device **[MD]** is made only for multiple uses on a single patient . If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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circaid® juxtafit®

Intended purpose

This compression product serves to compress the leg in patients with venous and lymphatic diseases .

Indications

- Lymphoedema
- Lipoedema
- Chronic venous insufficiency
 - Varicose veins
 - Lipodermatosclerosis
 - Venous ulcer
- Deep vein thrombosis/thromboprophylaxis
- Acute deep vein thrombosis
- Post-thrombotic syndrome
- Post sclerotherapy
- Dependent oedema

Contraindications

- Advanced peripheral arterial disease (if one of these parameters is present: ABPI < 0.5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot). When non-elastic materials are implemented, the patient can still try out compression clothing if they have ankle arterial pressure between 50 and 60 mmHg, under close clinical supervision.
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Untreated infection of the leg and/or foot
- Any circumstances under which an increase in venous or lymphatic reflux is undesirable.
- Suspected or known untreated acute deep vein thrombosis

Risks/side effects

- Slight or moderate peripheral arterial disease
- Decreased sensation – the patient must be able to feel the degree of compression
- Sensitivity to one of the used materials
- Use in children and patients dependant on assistance – patients must be able to feel and communicate the applied degree of compression

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Recommendations for wearing

Note: Unless your treating physician orders or recommends otherwise, we recommend wearing the legwear together with the respective footwear included in the set. If in doubt, please consult your treating physician.

Always make sure that the correct pressure is applied for compression. At night, slightly loosen the straps of the compression product. Every vertical BPS line on the straps can be at a different spot on the leg according to the size and shape of the leg. They do not have to be aligned to the front centre of the leg or to another defined area of the leg. In case of pain, reduce the degree of compression or remove the compression product. Notify your doctor. If applied correctly, the straps overlap slightly so that there are no gaps. Only tighten the straps in the intended places. The hook-and-loop fastener must never come into contact with skin!

Fitting instructions

- Put on the circaid undersleeve (understocking). The thick end of the understocking should cover the heel. Make sure there are no creases in the understocking. (Fig. 1)
- Position the compression product above the ankle and below the back of the knee. The black side of the compression product must face towards the leg. Unroll the two top straps and fasten them loosely on the opposite side to fix the compression product in the correct position. This makes adjusting the compression product easier. (Fig. 2)
- Unroll the two lowest straps and fasten the lowest strap on the opposite side of the compression product with a tight, but comfortable, degree of compression. (Fig. 3)
- Before fastening the second strap from the bottom, loosen the third strap from the bottom. Now fasten the second strap. Repeat this process with the remaining straps, working towards the knee. (Fig. 4)
- The compression product should fit snug and smoothly around the leg and the applied degree of compression should be tight but comfortable. Make sure that all straps have been fastened in an alternating order. (Fig. 5)
- Check the compression product for gaps and creases after all straps have been fastened. If necessary, adjust them. (Fig. 6)

Note: After adjusting the correct pressure using the BPS card, you can fold the lower and upper end of the circaid undersleeve (understocking) over the circaid compression product.

Putting on the circaid compressive sock

First fold the lower part of the understocking over the circaid compression product. Put on the compressive sock and make sure it overlaps with the product used on the leg. Check the compressive sock for creases and if necessary, smooth. If there is additional or exclusive use of another foot option, please pay attention to the notes in the respective instruction manual.

Note: Please roll up the straps when removing, the hook-and-loop fasteners must not be open.

Using the built-in-pressure-system (BPS)

The set compression value can vary due to movement, specific material characteristics and individual body anatomy.

- Pick up the BPS card.
- Identify the degree of compression prescribed for you on the colour coded card.
- Start with the lowest strap on the compression product. Place the black triangle on the BPS card against one of the BPS lines on the lowest strap.

- Check where the second BPS line of the strap meets the compression areas on the BPS card (A).
- If the BPS line is below (B) or above (C) the desired compression area, loosen and fasten the strap again until reaching the desired compression (A). If there is too much elasticity, completely unfasten the strap before tightening again.
- Repeat steps 3 - 5 with all straps along the leg. Make sure that all straps are adjusted to the same level of compression to achieve a decreasing rate of compression from distal to proximal. In order to maintain the range of compression throughout the day, if necessary check the distances between the BPS lines with the BPS-card. In the event of possible changes, please make adjustments back to the desired compression range.

Note: Please make sure to always cover the entire limb with the medical compression product.

Care instructions

This compressive product can be washed on a delicate cycle and dried in a tumble dryer set to a low heat. Please roll up the straps for laundry, the hook-and-loop-fasteners must not be open. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

-  Delicate wash cycle
-  Do not bleach
-  Tumble dry at low temperature
-  Do not iron
-  Do not dry clean
-  Do not wring

Storage instructions

Please ensure that the product is dry before storing it and store it out of direct sunlight.

Product life

6 months – juxtafit essentials

12 months – juxtafit premium

Due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product).

Material composition

Polyurethane, nylon, elastane

Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please notify your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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