**medi GmbH & Co. KG**  
**Medicusstraße 1**  
**95448 Bayreuth**  
**Germany**  
**T +49 921 912-0**  
**www.medi.de**

**medi Australia Pty Ltd**  
**Unit 4/13-21 Hallmark Street**  
**Pendle Hill NSW 2145**  
**Australia**  
**T +61 2 9890 8696**  
**F +61 2 9890 8439**  
**sales@medi.au**  
**www.mediaustralia.com.au**

**medi Austria GmbH**  
**Adamgasse 16/7**  
**6020 Innsbruck**  
**Austria**  
**T +43 512 57 95 15**  
**F +43 512 57 95 15 45**  
**vertrieb@medi-austria.at**  
**www.medi-austria.at**

**medi Belgium NV**  
**Staatsbaan 77/0099**  
**3945 Ham**  
**Belgium**  
**T: + 32-11 24 25 60**  
**F: +32-11 24 25 64**  
**info@medibelgium.be**  
**www.medibelgium.be**

**medi Canada Inc / médi Canada Inc**  
**104-1375 Lionel-Boulet,**  
**Varennnes, Québec,**  
**QC Canada J3X 1P7**  
**T +1 450-583-3317 / +1 800 361 3153**  
**F +1 888-583-6827**  
**service@medicanada.ca**  
**www.medicanada.ca**

**MAXIS a.s.,**  
**medi group company**  
**Slezská 2127/13**  
**120 00 Prague 2**  
**Czech Republic**  
**T: +420 571 633 510**  
**F: +420 571 616 271**  
**info@maxis-medica.com**  
**www.maxis-medica.com**

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### Important notes

This medical device **MD** is made only for multiple uses on a single patient. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If you experience excessive pain or any discomfort whilst wearing, remove the product and immediately contact your doctor or medical supply store. Only wear the product on unbroken skin or over injured skin that has been medically treated, not directly on injured or damaged skin, and only after previous consultation with a medical professional.

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## **circaid® juxtafit® essentials arm**

### Intended purpose

The compression system is designed to provide compression  to the arm for patients with venous and lymphatic disorders.

### Indications

- Lymphedema
- Lipedema
- Dependent oedema
- Acute shoulder-arm vein thrombosis

### Contraindications

- Uncontrolled Congestive Heart Failure
- Septic Phlebitis
- Any circumstance where increased venous and lymphatic return is undesirable
- Suspected or known untreated acute deep vein thrombosis

## Risks / Side Effects

- Slight or moderate peripheral arterial disease
- Impaired sensation – patient must be able to feel the degree of compression
- Sensitivity to one of the used materials (may cause skin irritation or itching)
- Paediatric and assisted use – patients must be able to sense and communicate applied pressure

## Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that they have received briefing on this from healthcare professionals.

Patient target group: Healthcare professionals should provide care to adults and children, applying the available information on the measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

## Note on wearing

When applied correctly, the straps overlap slightly so that there are no gaps. Only tighten the straps in the intended places. The hook-and-loop fastener must never come into contact with the skin. The compression product should feel firm, but comfortable.

## Fitting instructions

- You may wish to wear a circaid undersleeve. The circaid juxtafit essentials arm can also be worn directly against the skin.
- Slide the circaid juxtafit essentials arm onto the arm. The top of the armsleeve should stop comfortably below the axilla. The elastic material should be on the inner side of the arm with the elbow section of the sleeve centered with the elbow crease.
- Unroll the band closest to the wrist, pull it across the elastic section of the sleeve and attach it onto the opposite side with firm and comfortable compression. Repeat this process with every band continuing up the arm. The garment should lay flat and wrinkle-free against the arm. Adjust each band to a firm and comfortable compression level. Ensure that the bands are placed in an alternating order.

**Note:** Once all bands are secured inspect the garment for any gaps or creases. Adjust the bands as necessary. The circaid juxtafit essentials elbow strap will be pre-applied in its proper position for ease of

use however some adjusting to maximize comfort and coverage may be needed when wearing the arm garment.

Please note that the limb is always treated in its entirety with medical compression. This means that compression products for the arm begin at the hand or the wrist, and some extend up to just below the elbow or beyond to below the axilla.

## Using the built-in-pressure system (BPS)

The built-in-pressure system (BPS) is implemented on the circaid product with the BPS card and the BPS lines and is used to set, check and adjust the defined compression range of the circaid compression product. You can use the BPS card to set the prescribed compression range.

- Put on the circaid compression product so that it fits tightly, but comfortably (see fitting instructions).
- Pick up the BPS card
- Identify the degree of compression prescribed for you on the colour coded card
- Start with the lowest strap on the compression product. Place the black triangle on the BPS card against one of the BPS lines; the second BPS line is positioned on the prescribed degree of compression.
- If the BPS line is below (B) or above (C) the desired compression area, loosen and fasten the strap again until reaching the desired compression (A).
- Repeat this process for all straps from distal to proximal.

Always make sure that all straps are set to the same, prescribed level of compression to ensure graduated compression.

In order to maintain the prescribed compressive pressure throughout the day, check the distances between the BPS lines with the BPS card and adjust the straps at intervals throughout the day, if necessary. Always make sure that the correct, prescribed pressure is applied for compression.

Note: The BPS lines for individual straps can generally be found at different points; it is not necessary to align the BPS lines in one row lengthwise or in any other specific position.

**Note:** If using multiple circaid juxtafit made-to-measure garments, only lightly apply bands overlapping the legging below to ensure gradient compression. Set the remaining bands to firm and comfortable compression or prescribed compression level.

## Care instructions

This compressive product can be washed

on a delicate cycle and dried in a tumble dryer set to a low heat. Please roll up the straps for laundry, the hook-and-loop-fasteners must not be open. We recommend using a mesh laundry bag. To extend the life of the product, hand wash and drip dry.

-  Delicates wash cycle
-  Do not bleach
-  Tumble dry at low heat
-  Do not iron
-  Do not dry-clean
-  Do not wring out

### Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

### Lifespan

6 months – due to material wear and/or slackening, medical efficacy can only be guaranteed for a defined lifespan. This assumes correct handling of the product (e.g. during care and when putting on and taking off the product). Unless otherwise prescribed by your doctor, wear the appliance daily from morning to night.

### Material composition

Nylon, elastane, polyurethane

### Disposal

Dispose of properly after use.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR).

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