Maxis®


Návod k použití / Návod na použitie / Instructions for use /
Gebrauchsanweisung / Mode d'emploi / Instrucciones de uso /
Instruções de uso / Istruzioni d'uso / Gebruiksaanwijzing /
Brugsanvisning / Bruksanvisning / Bruksanvisning /
Käyttöohje / Naudojimo vadovas / Kasutusjuhend / Instrukcja
użytkowania / Használati utasítás / Navodila za uporabo /
Инструкция по применению / Kullanma talimatları / تعليمات الاستخدام

...we care about you


1

2

3

4

5

6

B

1

2

3

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ENGLISH

INSTRUCTIONS FOR USE

The instructions apply to the following versions:

MICRO, BRILLANT, COMFORT COTTON, PROMAX STANDARD, made-to-measure products

Dear user,
thank you for your trust. We are pleased that you have opted for our medical compression stockings. To ensure the effectiveness of your compression therapy, we ask that you take some time to read these instructions carefully. Keep them in a safe place so that you can find important information later if needed. If you have any questions, please contact your doctor or healthcare staff you trust.

1. Purpose – what are medical compression stockings used for?

Medical compression stockings are an integral part of compression therapy. Compression therapy is part of the treatment for venous, lymphatic and mixed edema and is one of the irreplaceable methods of conservative care.

2. Performance characteristics – how do health compression stockings work?

The mechanism of action is to reduce capillary filtration into the interstitium, assist in the resorption of tissue fluid (and crystalloids) into the venous capillary arm, accelerate lymph flow through the lymphatic system, and enhance the venous-lymphatic-muscle pump effect. This spectrum of effects minimizes or completely eliminates the adverse effects of venous reflux, venostasis and lymphostasis.

3. Indications – what symptoms are medical compression stockings used for?

After diagnosis you have received our medical compression device.
For optimal treatment results, follow the instructions of your doctor or healthcare staff. Please note that each stocking has different indications and areas of use due to the nature of the material and the available compression classes. Your medical compression stockings can generally be used for the following indications:

Chronic venous diseases

- improvement of venous symptoms
- improvement of the quality of life in chronic venous diseases
- prevention and therapy of venous edema
- prevention and therapy of venous skin changes
- eczema and pigmentation
- dermatoliposclerosis and white atrophy
- treatment of leg ulcers
- treatment of mixed (arterial and venous) ulcer cruris

(with regard to contraindications: see point 4)

- prevention of venous leg recurrence
- reduction of leg ulcer pain
- varicose veins
- initial phase after varicose veins treatment
- functional venous insufficiency (for obesity, sedentary and standing jobs)
- venous malformations

Thromboembolic venous diseases

- superficial venous thrombosis
- deep vein thrombosis
- thrombosis
- post-thrombotic syndrome
- prophylaxis of thrombosis in mobile patients

Edema

- lymphedema
- edema in pregnancy
- post-traumatic edema
- postoperative edema
- postoperative reperfusion edema
- cyclic idiopathic edema
- lipedema
- congestion due to immobility (arthrogenic congestion syndrome, paresis and partial limb paresis)
- working edema (sedentary and standing jobs)
- drug-induced edema

Other indications

- obesity with functional venous insufficiency
- inflammatory dermatosis of the feet
- nausea, dizziness during pregnancy
- problem with overburdening during pregnancy
- post-burn condition
- treatment of scars

Please note the following information:

Patient-specific factors such as body weight, type and seve-

ity of swelling, and the nature of connective tissue also play a role in the selection of medical compression stockings.

Therefore, our recommendations:

- The higher your body weight,
 - The greater the tendency of swelling,
 - The more severe the disease,
 - The softer the connective tissue,
- the thicker the stocking material should be!

Your doctor or medical staff you trust will help and advise you on the correct choice of medical compression stockings. This will ensure that your individual needs are taken into account – for greater wellbeing and therapeutic success.

4. Contraindications – when should medical compression stockings not be worn?

Medical compression stockings must not be worn in the following cases:

- advanced peripheral arterial obliteration (if one of the following applies: ABPI < 0.5, ankle arterial pressure < 60 mm Hg, finger pressure < 30 mm Hg or TcPO2 < 20 mm Hg instep).
- Using inelastic materials, a dose of compression at ankle arterial pressure between 50 and 60 mm Hg can be tested under clinical control.
- decompensated cardiac insufficiency (NYHA III + IV)
- septic phlebitis
- phlegmasia coerulea dolens

In the following cases, the therapeutic decision should be carefully considered (in view of the benefits and risks and the choice of the most appropriate compression):

- severe wet dermatosis
- intolerance to compression material
- severe limb sensitivity disorders
- advanced peripheral neuropathy (e.g. in diabetes mellitus)
- primarily chronic polyarthritis

If you are not sure if any of these apply to you, talk to your doctor or a healthcare professional you trust.

The following risks and side effects may occur:

In exceptional cases, medical compression stockings can cause, especially in the case of improper handling

- skin necrosis and
 - pressure damage of peripheral nerves.
- Sensitive skin may be prone to itchiness, scaling and signs of inflammation. Therefore, it is reasonable to take care of your skin during compression care. Please follow our important guidelines (points 7 and 8).

The following symptoms should result in the immediate removal of the compression sleeve and consultation with the attending physician or a healthcare professional:

- blue or white toes, discomfort and numbness, increasing pain, shortness of breath and sweating, acute movement restriction

5. Intended target user group

The intended users include doctors, healthcare professionals, nursing staff ensuring correct application of the medical compression sleeve to patients, and the patients themselves. The correct choice of medical compression sleeve (measurement, size determination, function, indication in adults and children, compression class) is their responsibility and should be based on available information from the manufacturer.

6. Wearing and use time – how long can I wear compression stockings?

Unless your doctor specifies otherwise, wear your compression stockings from morning to evening for an optimal therapeutic effect. If in exceptional cases (e.g. a long-haul flight) the wearing time is significantly exceeded, you have the possibility to slightly adjust the compression stockings several times.

For the stockings to work optimally, a precisely graduated pressure curve is a prerequisite. Everyday wearing and washing can reduce the necessary therapeutic pressure and elasticity of your compression stockings over time. The recommended usability is therefore 6 months. To obtain the best pressure, follow our guidelines and recommendations (points 8 and 10). Do not hesitate to contact us. We will inform you about the possibilities of ensuring effective and long-term treatment of diseases of the venous or lymphatic system using medical compression stockings.

7. Instructions for use – what should I know?

- Medical compression stockings should only be distributed by trained personnel and accompanied by instructions for use and care (see points 8 and 10).
- Should you experience pain or excessive skin irritation during wear, remove your stockings immediately and consult your doctor or a healthcare professional.
- Appropriate skin care is essential for compression therapy. In particular, creams or ointments containing fat, but also soap residues, can cause skin irritation and material wear, and thus reduce the effectiveness of medical compression stockings. That's why we offer complementary skin care products that are specially adapted for compression stockings. Ask for advice at specialized points of sale.
- In silicone banded compression stockings, skin irritation may occur in people with sensitive skin (due to a combination of sweat and mechanical stress). To avoid this, it is recommended to slightly move the band several times during the day (see point 6). Slipping of thigh stockings can also be caused by a large amount of hairs. In similar cases, we recommend the use of a special adhesive for a better grip of the medical compression stockings. If necessary, ask for advice at specialized points of sale.
- Sharp finger and toe nails, defective shoes, jewelry, as well as rough and calloused heels can cause damage to the knit.
- Never repair stockings by yourself or have them repaired by another service provider, as this reduces the quality, safety and efficiency of the medical device. Information on repair, if possible, will be provided to you by our salespeople or our customer service.
- In case of complaints about products such as knit damage, please contact your specialized shops. Only serious events that may lead to a significant deterioration of health or death shall be reported to the manufacturer and the competent state authority. Major events are defined in Article 2 no 65. of Regulation (EU) 2017/745 (MDR).
- A textile label sewn onto the stockings plays an important part in identifying the product. Under no circumstances should you remove this label from your medical compression stocking, cutting it will void your warranty, repair, or replacement.

8. Donning instructions – what should I be careful about when putting on my medical compression stockings?

Always don the product on a dry leg, preferably in the morning before putting pressure on the limb.
If you take a shower in the morning and then take care of your skin (see also point 10 – Instructions for care), wear a short time until your skin dries properly. This makes donning much easier. Also prepare your textile or rubber gloves, which are available at specialized points of sale. You can use them later to make it easier to spread the stocking around your leg or to get started at the beginning. The gloves fit particularly well and make donning easier.
When putting the stockings on, make sure your nails are rounded and smooth. Well-kept heels and toenails are also important to prevent damage to the knit. Put on your jewelry or wristwatch only after putting on your stockings. To put your medical compression stockings on comfortably, sit down and place your foot on an elevated surface, such as a small stool.

Donning open-toed medical compression stockings
The procedure for putting on open-toed products is shown in the front part of the instructions (fig. A).

- For easier donning of open-toed products, we have included a smooth donning aid which goes on your toes and heel.

CZ
SK
EN
DE
FR
ES
PO
IT
NE
DA
SWE
NO
FI
LT
ET
PL
HU
SL
RU
TU
AR

2. Put the stocking on inside out up to the ankle area. Insert your thumbs inside the stocking and stretch the stocking over your toes.
3. Pull the stocking over the heel. The heel is the area where you need to develop more force to pull up the stocking.
4. The highest compression pressure is in the ankle area, the product is narrowest here. Gradually pull up the next part of the stocking by scooping up the knit.
5. Proceed further and pull the stocking towards the knee, or higher, depending on the type of stocking. Finally, smooth out any creases and unevenness in the direction from the ankle with your palms until you feel the tight compression effect of the product over the entire surface of your leg.
6. Finally, pull the smooth donning aid out through the open toe of the stockings.
- Initially, you may find donning difficult, but with repeated use you will gain the skill you need and the process will become easier for you.

Donning closed-toe medical compression stockings

The donning procedure is illustrated in the front part of the instructions (fig. B).

1. Don products with a closed toe by inserting your thumbs inside the stocking and pulling the stocking gradually over your toes and over the heel. The heel is an area where you need to develop more force to pull up the stockings.
2. The highest compression pressure is in the ankle area, the product is narrowest here. Gradually pull up the next part of the stocking by scooping up the knit.
3. Proceed further and pull the stocking towards the knee, or higher, depending on the type of stocking. Finally, smooth out any creases and unevenness in the direction from the ankle with your palms until you feel the tight compression effect of the product over the entire surface of your leg.

The threading procedure is illustrated in the front part of the instructions for use (fig. B: Donning closed-toe medical compression stockings).

Tip: Ask for a practical metal donning aid – a medical compression stocking donning aid.

9. Storage guidelines – how and how long can I store my medical compression stockings?

The products should be stored in a dry place and protected from direct sunlight and heat.

-   The hourglass symbol on the box label indicates the latest usability date for the product.
-  This is the maximum shelf life of serial products in the length of 42 months. This period already includes a maximum practical usability of 6 months (see point 6).
- Made-to-measure products are designed for immediate use.

10. Care guidelines – what do I have to look out for when caring for my medical compression stockings?

Wash your stockings every day after wearing. We recommend handwashing detergent for gentle care. Alternatively, use a mild laundry detergent (delicate wash at 40° C) without an optical brightener and softener.

Tip:

- We recommend turning your compression stockings inside out when washing. Wash them separately or with clothes of the same color.
- A washing bag better protects the stocking material.
- Do not dry clean medical compression stockings or wash them using household detergents.
- After washing, you can roll the stockings into a towel and squeeze. Do not leave wet stockings in a wet towel, air dry them. Under no circumstances should they be dried in the sun or on a radiator. Alternatively, you can delicately tumble dry your medical compression stockings.

**Wash up to 40° C (delicates)**

**Do not bleach**

**Do not iron**

**Do not dry clean**

**Low tumble dry temperatures (delicates)**

11. Material composition – what are your compression stockings composed of?

Compression stockings MICRO, BRILLANT, COTTON and PROMAX STANDARD contain polyamide and elastane. COTTON stockings also contain cotton. More detailed information can be found on the textile label sewn on your compression stockings.

 Without latex: There is no rubber from the milk juice of tropical plants processed in this product.

- MAXIS MICRO Compression Stockings include MICRO PROTECTION. It is an antibacterial and antifungal treatment that prevents the growth of fungi, bacteria and other undesirable microorganisms.
- MAXIS BRILLANT compression stockings incorporate a special High IQ® COOL COMFORT treatment for better moisture management.
- MAXIS COTTON compression stockings contain microcapsules with an antibacterial effect in combination with Aloe Vera extract.
- PROMAX STANDARD compression stockings include SOFT PROTECTION. This is an antibacterial finish that makes the fine and thin stockings fresh and breathable. It also contributes to better moisture removal from the body.

12. Disposal – how do I dispose of my medical compression stockings?

Please dispose of your medical compression stockings in your normal household waste. There are no special criteria for disposal.

13. Explanation of symbols

**Medical Device**

**Unique Device Identifier**

**Single patient – multiple use (see also point 6: Wearing time and usability)**

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