
mediven[®] thrombexin[®] 18

mediven[®] thrombexin[®] 18 Anti-slip

Gebrauchsanweisung. Instructions for use. Инструкции за употреба. Brugsvejledning. Kasutusjuhend. Käyttöohje. Mode d'emploi. Οδηγίες χρήσης. 使用説明書. Istruzioni per l'uso. Uputa za uporabu. Lietošanas instrukcija. Naudojimo instrukcija. 使用说明. Gebruiksaanwijzing. Bruksanvisning. Instrukcja załadania. Manual de instrucções. Instrucţiuni de utilizare. Инструкция по использованию. Bruksanvisning. Návod na použitie. Navodila za uporabo. Instrucciones de uso. Návod k používání. Kallanma kılavuzu. تعليمات الاستعمال. Інструкція з використання. Használati útmutató. הוראות שימוש. Arahan penggunaan. Uputstvo za upotrebu.


1


Medizinprodukt / Medical product / Медицинско изделие / Medicinsk produkt / Meditsiiniseade / Lääkinnällinen tuote / Produit médical / Ιατροτεχνολογικό προϊόν / 醫療產品 / Presidio medicale / Medicinski proizvod / Mediciniska ierice / Medicinos priemonė / 医疗器械 / Medisch product / Medicinsk produkt / Produkt medyczny / Dispositivo médico / Produs medical / Медицинское Изделие / Medicinteknisk produkt / Zdravotnícky produkt / Medicinski pripomoček / Producto médico / Zdravotnícký prostředek / Tibbi ürün / Медицинский виріб / منتج طبي / Orvos-technikai eszköz / מיכשור רפואי / Medicinski proizvod


Anzahl der Wiederverwendungen minus eins / Number of times product can be reused minus one / Броу на повторните употреби минус едно / Antal gange gentaget brug minus en / Korduvkasutuste arv miinus üks / Uudelleenkäyttöjen lukumäärä miinus yksi / Nombre de réutilisations moins un / Αριθμός εκ νέου χρήσεων μείον ένα / 重複使用次數減一 / Numero di riutilizzi meno uno / Broj ponovnih korištenja minus jedan / Atkārtoto izmantojumu skaits minus viens / Panaudojimo kartų skaičius atėmus vieną / 复用次數減1 / Aantal keren dat een product kan worden hergebruikt min één / Antall gjenbruk minus ett / Liczba określająca, ile razy można założyć pończochę, pomniejszona o jeden / Número de vezes que o produto pode ser reutilizado menos um / Număr de reutilizări minus una / Количество повторных циклов „использование-стирка“ минус один / Antal gånger som produkten kan återanvändas minus en / Počet opätovných použití minus jedna / Število ponovnih uporab, zmanjšano za ena / Número de repeticiones de uso menos uno / Počet opakovaných použití minus jedna / Yeniden kullanım sayısı eksi bir / Кількість разів повторного використання продукту, мінус один / عدد مرات إعادة الاستخدام ناقص واحد / מספר השימושים החוזרים פחות אחד / Broj ponovnih upotreba minus jedan


2


3


4


5

6

7

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Intended purpose

mediven® thrombexin® 18:

Round knitted clinical compression  stocking for compression of the lower extremities, mainly for prevention of diseases of the venous system.

Product features

These products work using a graduated pressure pattern (decreasing from the end of your limb towards your heart) to exert compression on your limbs.

Indications

mediven® thrombexin® 18:

- Physical thrombosis prophylaxis in immobile patients - pre-, intra- and post-operative
- Prevention and therapy of post-operative and post-traumatic oedema

Contraindications

Advanced peripheral arterial occlusive disease (if one of these parameters is present: ABPI < 0,5, ankle arterial pressure < 60 mmHg, toe pressure < 30 mmHg or TcPO₂ < 20 mmHg on dorsum of foot), Decompensated heart failure (NYHA III + IV), Septic thrombophlebitis, Phlegmasia cerulea dolens, Massive leg oedema

Take especial care if the patient has

Severe weeping dermatosis, Intolerance to compressive material, Severe paraesthesia in the limbs, Advanced peripheral neuropathy (e.g. in diabetes mellitus), Primary chronic polyarthritis

Risks and side effects

The products may incur the following, especially if handled incorrectly: skin necrosis and pressure damage to peripheral nerves. If your skin is sensitive, the compressive equipment may cause itchiness, peeling or signs of inflammation. Adequate skincare under the compression stocking is therefore vital. Please comply with the instructions for use and inspect your skin each day, checking that your circulation is optimal. If you experience the following symptoms, the compression stocking should be removed immediately and clinical advice should be sought: Change in toe colour to blue or white, new paraesthesia or numbness, an increase in pain, shortness of breath, sudden sweating or acute limitation of mobility.

Intended users and patient target groups

Healthcare professionals and patients, including persons who play a supporting role in care, are the intended users for this device, provided that the persons supporting carer has received a briefing on this from healthcare professionals. Healthcare professionals should provide care to adults and children, applying the available information on measurements/sizes and required functions/indications, in line with the information provided by the manufacturer, and acting under their own responsibility.

Instructions for use

Optimal care can only be provided if the measurements are taken and applied correctly.

mediven® thrombexin® 18:

The circumference of the ankle (CB) dictates the size that must be selected. For the knee stocking, the measurement of the calf circumference (CC) confirms the correct size. For thigh stockings, the selection of the size is confirmed by the measurement of the thigh circumference (CG), or this may indicate that an extra-wide stocking is needed. In the interests of better handling and legibility, measurements for garments are given in whole centimetres. Any potential gaps that may result from this

are covered by measurement and production tolerances. The internationally recognised rules for rounding apply for the measurements (round down for numbers up to 4/10, round up from 5/10). This was taken into account when designing the paper measuring tapes provided by the manufacturer.

The products may only be used on intact skin or if an appropriate wound dressing is in place, and are suitable for use by multiple patients, provided that they are cleaned in advance. Contact with chemical substances (e.g. benzenes, methanols, salicylic acids in high concentrations) may impair the compressive function of the products. Oils and ointments should not be applied to the limb when compression stockings are used. Compression stockings with a silicone top band can cause skin irritation in people with sensitive skin (caused by perspiration in combination with mechanical stress). In order to prevent this, it is recommended that the top band is repositioned slightly several times throughout the day. If you experience pain or increased skin irritation, take off the compression stocking immediately and consult your doctor or carer.

Length of time that compression clothing should be worn and used for

Unless otherwise stipulated by your doctor, you can wear the compression stocking day and night in order to elicit optimal therapeutic results. The compression stockings should be changed as often as necessary. Under normal circumstances, this should be every 3 days, at the least. The compression stockings can be laundered a maximum of 10 times.

Instructions for care

Launder according to the requirements of the relevant Medical Device Directive and the relevant national requirements for laundering medical devices by regarding the following preconditions:

- Wash at temperatures of up to 95°C
- Do not bleach
- Dry at temperatures of up to 80°C
- Do not iron
- Do not dry-clean

Damaged compression stockings must be removed from circulation and can no longer be used on patients for the set indications.

Storage

The products must be stored in a dry place and out of direct sunlight.

Disposal

Dispose of properly after use.

Complaints

In the event of any complaints regarding the product, such as damage to the fabric or a fault in the fit, please contact the manufacturer. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer and the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of Regulation (EU) 2017/745 (MDR).

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medi UK Ltd.
Cosanum AG